# **ISOBUTYL ACETATE**



## 1. PERFORMANCE

1) Measuring range 0.01-1.4 % Number of pump strokes  $2(200 \text{m} \ell)$ 

2) Sampling time : 4 minutes/2 pump strokes

3) Detectable limit  $\therefore$  10 ppm 4) Shelf life  $\therefore$  3 years 5) Operating temperature  $\therefore$  0  $\sim$  40 °C

6) Reading : Graduations printed on the tube are calibrated by Methyl ethyl ketone at 2 pump strokes and Isobutyl acetate concentration is determined by using a conversion chart.

7) Colour change : Orange → Pale brownish green

## 2. RELATIVE STANDARD DEVIATION

RSD-low: 15% RSD-mid.: 10% RSD-high: 5%

## 3. CHEMICAL REACTION

Potassium dichromate is reduced CH<sub>3</sub>CO<sub>2</sub>CH<sub>2</sub>CH (CH<sub>3</sub>)<sub>2</sub> + Cr<sup>6+</sup> + H<sub>2</sub>SO<sub>4</sub>→ Cr<sup>3+</sup>

### 4. CALIBRATION OF THE TUBE

GAS CHROMATOGRAPHY

#### 5. INTERFERENCE AND CROSS SENSITIVITY

Substance	Interference	ppm	Coexistence
Acetylene		3%	Whole reagent is changed to Pale brown.
Propane		0.2%	"
Other organic gases or vapours except Halogenated hydrocarbons	Similar stain is produced.	50	Higher readings are given.





No.139SB Tube reading (%)