
ISOPENTYL ALCOHOL



1. PERFORMANCE

1) Measuring range 5-100 ppmNumber of pump strokes $3(300 \text{m} \ell)$

2) Sampling time : 4.5 minutes/3 pump strokes

3) Detectable limit \therefore 2 ppm 4) Shelf life \therefore 2 years 5) Operating temperature \therefore 0 \sim 40 $^{\circ}$ C

6) Temperature compensation : Necessary (See "TEMPERATURE CORRECTION TABLE") 7) Reading : Direct reading from the scale calibrated by 3 pump strokes

8) Colour change : Yellow→Pale blue

2. RELATIVE STANDARD DEVIATION

RSD-low: 10% RSD-mid.: 10% RSD-high: 5%

3. CHEMICAL REACTION


Chromium oxide is reduced. (CH₃)₂CHCH₂OH + Cr⁶⁺ + H₂SO₄→Cr³⁺

4. CALIBRATION OF THE TUBE

GAS CHROMATOGRAPHY

5. INTERFERENCE AND CROSS SENSITIVITY

0 0 0 0										
Substance		ppm Interference		Coexistence						
Alcohols	FIG.1		Similar stain is produced.	Higher readings are given.						
Toluene	FIG.2	200	Whole reagent is changed to Pale blue.	"						
Hexane		1,000	The accuracy of readings is not affected.	The accuracy of readings is not affected.						
Trichloroethylene		"	"	"						
Ethyl acetate		"	"	"						



TEMPERATURE CORRECTION TABLE

Scale	True Concentration (ppm)						
Readings (ppm)	0 °C (32 °F)	10℃ (50°F)	20°C (68°F)	30°C (86°F)	40°C (104°F)		
100	_	_	100	80	70		
80	_	_	80	65	55		
60	_	83	60	50	45		
40	90	50	40	33	30		
20	32	23	20	16	15		
10	14	11	10	9	8		
5	5	5	5	5	5		